


L-Homoserine lactone

Code No.: **BIA-H2224**

Pack sizes: **1 mg, 5 mg**



Synonyms : L-Homoserine γ -lactone, L- α -Aminobutyrolactone, (3S)-3-Aminodihydro-2(3H)-furanone

Specifications

CAS #	: 2185-02-6
Molecular Formula	: C₄H₇NO₂
Molecular Weight	: 101.1
Source	: Synthetic
Appearance	: White solid
Purity	: >95% by HPLC
Long Term Storage	: -20°C
Solubility	: Soluble in methanol or DMSO

Application Notes

References

1.