


## Lagosin

Code No.: **BIA-L1613**

Pack sizes: **1 mg, 5 mg**



Synonyms : Fungichromin, 14-Hydroxyfilipin III, Cogomycin, Moldcidin B, A246, Pentamycin

## Specifications

CAS #	: 6834-98-6
Molecular Formula	: C <sub>35</sub> H <sub>58</sub> O <sub>12</sub>
Molecular Weight	: 670.8
Source	: <i>Streptomyces</i> sp.
Appearance	: Pale yellow to tan solid
Purity	: >95% by HPLC
Long Term Storage	: -20°C
Solubility	: Soluble in ethanol, methanol, DMF or DMSO.

## Application Notes

Lagosin is a pentaene antifungal produced by *Streptomyces*, first isolated in 1958 by researchers at MIT in the USA. The discovery was soon followed by several independent isolations as lagosin and cogomyin. Initially these metabolites were thought to be isomeric, but Pandey and colleagues at NCI definitively demonstrated they were identical. Structurally, lagosin is 14-hydroxyfilipin III and the most polar member of the filipin family of fungicides. Lagosin exhibits broad spectrum antifungal and antitumor activity and, like filipin, acts via interaction with cell membrane sterols.

## References

1. Fungichromin: Determination of the structure of the pentaene chromophore. Cope A.C. & Johnson H.E. J. Am. Chem. Soc. 1958, 80, 1504.
2. Physicochemical and biological comparison of polyene macrolide antibiotics fungichromin, lagosin and cogomycin. Pandey R.C. et al., J. Antibiot. 1982, 35, 988.