


Pseudoerythromycin A enol ether

Code No.: BIA-P1349

Pack sizes.: 1mg, 5mg



Synonyms:

-

Specifications

CAS #	: 105882-69-7
Molecular Formula	: C37H65NO12
Molecular Weight	: 715.9
Source	: -
Appearance	: White solid
Purity	: >98% by HPLC
Long Term Storage	: -20°C
Solubility	: Soluble in ethanol, methanol, DMF or DMSO. Good water solubility.

Application Notes

Pseudoerythromycin A enol ether is a degradation product of erythromycin formed by a complex internal rearrangement of erythromycin A on exposure to neutral to weakly alkaline conditions. The C6-OH forms an internal enol ether with the C9 ketone of erythromycin, while the C11-OH attacks the carbonyl of the lactone to reduce the macrocycle from a 14- to an 11-membered macrolide. Synthetically, pseudoerythromycin A enol ether is prepared by reacting erythromycin enol ether with carbonate. Pseudoerythromycin A enol ether is devoid of antibiotic activity but is an important analytical standard for erythromycin A stability studies.

References

1. Study of the stability of erythromycin in neutral and alkaline solutions by liquid chromatography on poly(styrene-divinylbenzene). Paesen J. et al., Int. J. Pharm. 1994, 113, 215.
2. Synthesis of ring-contracted derivatives of erythromycin. Kirst H.A. et al., J. Org. Chem. 1987, 52, 4359.