


Psychotridine

Code No.: BIA-P1960

Pack sizes.: 1mg, 5mg



Synonyms:

-

Specifications

CAS #	: 52617-25-1
Molecular Formula	: C55H62N10
Molecular Weight	: 863.15
Source	: -
Appearance	: Off-white solid
Purity	: >95% by HPLC
Long Term Storage	: -20°C
Solubility	: Soluble in methanol or DMSO

Application Notes

Analgesic, NMDA receptor antagonist. Hypotensive; histamine blocking activity; no analgesic, antipyretic, diuretic or anticonvulsant activity

References