


## L-Phe-L-Pro

Code No.: BIA-P2369

Pack sizes.: 5mg, 25mg



### Synonyms:

L-Phenylalanyl-L-proline

## Specifications

CAS #	: 7669-65-0
Molecular Formula	: C <sub>14</sub> H <sub>18</sub> N <sub>2</sub> O <sub>3</sub>
Molecular Weight	: 262.3
Source	: -
Appearance	: <b>White solid</b>
Purity	: <b>&gt;95% by HPLC</b>
Long Term Storage	: <b>-20°C</b>
Solubility	: <b>Soluble in methanol or DMSO</b>

## Application Notes

Amino acid

## References