


L-Phe-L-Val

Code No.: BIA-P2370

Pack sizes.: 5mg, 25mg



Synonyms:

L-Phenylalanyl-L-valine, Phenylalanylvaline

Specifications

CAS #	: 3918-90-9
Molecular Formula	: C₁₄H₂₀N₂O₃
Molecular Weight	: 264.32
Source	: -
Appearance	: White solid
Purity	: >95% by HPLC
Long Term Storage	: -20°C
Solubility	: Soluble in methanol or DMSO

Application Notes

Amino acid

References