


Code No.: BIA-P2467

Pepstatin A

Pack sizes.: 0.1mg, 0.5mg



Synonyms:

Ahpatinin C, NSC 272671, Pepsin inhibitor S 735A Pepstatin, Procidin S 735A

Specifications

CAS #	: 26305-03-3
Molecular Formula	: C34H63N5O9
Molecular Weight	: 685.89
Source	: -
Appearance	: White solid
Purity	: >95% by HPLC
Long Term Storage	: -20°C
Solubility	: Soluble in methanol or DMSO

Application Notes

-

References