


Setosusin

Code No.: BIA-S1084

Pack sizes.: 0.5mg, 2.5mg



Synonyms:

-

Specifications

CAS #	: 182926-45-0
Molecular Formula	: C₂₉H₃₈O₈
Molecular Weight	: 514.6
Source	: -
Appearance	: White Lyophilisate
Purity	: >99% by HPLC
Long Term Storage	: -20°C
Solubility	: Soluble in ethanol, methanol, DMF or DMSO. Poor water solubility.

Application Notes

Setosusin is an unusual bis-spiro metabolite originally isolated from *Corynascus setosus*. There is little published information on the biological profile of this metabolite.

References

1. Isolation of new tremorgenic metabolites from an ascomycete, *Corynascus setosus*. Fujimoto H. et al., Chem. Pharm. Bull. 1996, 44, 1843.