


Stauprimide

Code No.: BIA-S1226

Pack sizes.: 0.5mg, 2.5mg



Synonyms:
N-Benzoyl-7-oxostaurosporine

Specifications

CAS #	: 154589-96-5
Molecular Formula	: C₃₅H₂₈N₄O₅
Molecular Weight	: 584.6
Source	: -
Appearance	: Yellow orange solid
Purity	: >95% by HPLC
Long Term Storage	: -20°C
Solubility	

: Soluble in ethanol, methanol, DMF or DMSO. Limited water solubility.

Application Notes

Stauprimide is a semi-synthetic analogue of the staurosporine family of indolocarbazoles. Stauprimide was first published in 1994 as part of an extensive structure-activity investigation to improve the selective inhibition of protein kinase C as potential antitumor agents. Stauprimide increases the efficiency of the directed differentiation of mouse and human embryonic stem cells in synergy with defined extracellular signalling cues. Stauprimide interacts with NME2 (PUF) transcription factor to down-regulate c-Myc expression, leading to differentiation of stem cells.

References

1. Inhibitory activity and selectivity of staurosporine derivatives towards protein kinase C. Caravatti G. et al., *Bioorg. Med. Chem. Letters* 1994, 4, 399.
2. A small molecule primes embryonic stem cells for differentiation. Zhu S. et al., *Cell Stem Cell* 2009, 4, 416