


## Securinine

Code No.: BIA-S1947

Pack sizes.: 1mg, 5mg



### Synonyms:

Securinan-11-one, (-)-Securinine, NSC 107413, Securinin, l-Securinine

## Specifications

CAS #	: <b>5610-40-2</b>
Molecular Formula	: <b>C<sub>13</sub>H<sub>15</sub>NO<sub>2</sub></b>
Molecular Weight	: <b>217.26</b>
Source	: -
Appearance	: <b>Tan solid</b>
Purity	: <b>&gt;95% by HPLC</b>
Long Term Storage	: <b>-20°C</b>
Solubility	: <b>Soluble in methanol or DMSO</b>

## Application Notes

GABA antagonist, Antitumor

## References