


## Sartorypyrone D

Code No.: BIA-S1984

Pack sizes.: 0.1mg, 0.5mg



### Synonyms:

20-O-Deacetylsartorypyrone A

## Specifications

CAS #	: <b>1801167-48-5</b>
Molecular Formula	: <b>C<sub>26</sub>H<sub>38</sub>O<sub>4</sub></b>
Molecular Weight	: <b>414.58</b>
Source	: -
Appearance	: <b>White solid</b>
Purity	: <b>&gt;95% by HPLC</b>
Long Term Storage	: <b>-20°C</b>
Solubility	: <b>Soluble in methanol or DMSO</b>

## Application Notes

NADH-fumarate reductase inhibitor

## References