


Simvastatin acid

Code No.: BIA-S2271

Pack sizes.: 1mg, 5mg



Synonyms:

-

Specifications

CAS #	: 121009-77-6
Molecular Formula	: C ₂₅ H ₄₀ O ₆
Molecular Weight	: 436.6
Source	: -
Appearance	: White solid
Purity	: >95% by HPLC
Long Term Storage	: -20°C
Solubility	: Soluble in methanol or DMSO

Application Notes

Reversible inhibitor of HMGCR (HMG-CoA reductase), Cholesterol lowering

References