


Syringic acid

Code No.: BIA-S2355

Pack sizes.: 5mg, 25mg



Synonyms:

Cedar acid, Gallic acid 3,5-dimethyl ether, NSC 2129

Specifications

CAS #	: 530-57-4
Molecular Formula	: C9H10O5
Molecular Weight	: 198.17
Source	: -
Appearance	: White solid
Purity	: >95% by HPLC
Long Term Storage	: -20°C
Solubility	: Soluble in methanol or DMSO

Application Notes

-

References