


Shikimic acid

Code No.: BIA-S2407

Pack sizes.: 5mg, 25mg



Synonyms:

(-)-Shikimic acid, L-Shikimic acid, NSC 59257, Shikimic acid

Specifications

CAS #	: 138-59-0
Molecular Formula	: C7H10O5
Molecular Weight	: 174.15
Source	: -
Appearance	: White solid
Purity	: >95% by HPLC
Long Term Storage	: -20°C
Solubility	: Soluble in methanol or DMSO

Application Notes

-

References