


Sorokinianin

Code No.: BIA-S2565

Pack sizes.: 1mg, 5mg



Synonyms:  
(+)-Sorokinianin

## Specifications

CAS #	: <b>162616-73-1</b>
Molecular Formula	: <b>C<sub>18</sub>H<sub>28</sub>O<sub>4</sub></b>
Molecular Weight	: <b>308.41</b>
Source	: -
Appearance	: <b>Off-white solid</b>
Purity	: <b>&gt;95% by HPLC</b>
Long Term Storage	: <b>-20°C</b>
Solubility	: <b>Soluble in methanol or DMSO</b>

## Application Notes

-

## References