


Tetrahydroxyquinone

Code No.: BIA-T2306

Pack sizes.: 25mg, 100mg



Synonyms:

HPEK-1, Kelox, NSC 112931, THQ, Terasin

Specifications

CAS #	: 319-89-1
Molecular Formula	: C₆H₄O₆
Molecular Weight	: 172.09
Source	: -
Appearance	: Black solid
Purity	: >95% by HPLC
Long Term Storage	: -20°C
Solubility	: Soluble in methanol or DMSO

Application Notes

-

References