


L-Val-L-Ser

Code No.: BIA-V2367

Pack sizes.: 5mg, 25mg



Synonyms:

L-Valyl-L-serine, N-L-Valyl-L-serine, NSC 339921, Valylserine

## Specifications

CAS #	: <b>13588-94-8</b>
Molecular Formula	: <b>C<sub>8</sub>H<sub>16</sub>N<sub>2</sub>O<sub>4</sub></b>
Molecular Weight	: <b>204.22</b>
Source	: -
Appearance	: <b>White solid</b>
Purity	: <b>&gt;95% by HPLC</b>
Long Term Storage	: <b>-20°C</b>
Solubility	: <b>Soluble in methanol or DMSO</b>

## Application Notes

Amino acid

## References